ProphyCair & ProphyCair+

CE

9200-0727
version: 01

ENGLISH

User manual and technical description

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E-mail: info@linet.cz

<http://www.linet.cz>

ProphyCair & ProphyCair+
Mattress replacement system

Author: LINET, s.r.o.

Related links: www.linet.cz

Version: 01

Print date: 04/2016

Copyright © LINET spol. s r. o., 2016

Translation © LINET spol. s r. o., 2016

All rights reserved. All trademarks and brands are the property of the appropriate owners. The manufacturer reserves the right to changes in the contents of this manual that relate to the product's technical regulations. It is for this reason that the contents of this manual may indicate differences from the current manufacture of the product.

Content

1	Symbols	3
	Symbols and Labels on the Product.....	3
	Other Symbols.....	3
	Warning Notices	4
2	Contraindications.....	5
3	Technical Specifications.....	5
	Electrical Specifications	5
	Performance Specifications	5
	Mechanical Specifications	5
	Environmental Specifications.....	6
4	Safety Instructions	6
5	Standard and Regulations.....	7
6	Scope of Delivery	7
	Scope of Delivery.....	7
	List of Parts.....	7
7	Use.....	7
	Preparing the Bed for the Patient	7
	Putting the Patient into the Bed	7
8	Product Description - Control Unit.....	8
	Controls and Indicators.....	9
	Operating Instructions.....	9
9	Product Description - Mattress	11
	Mattress Weight Capacity.....	11
10	Unpacking the System.....	11
	Components Supplied	11
	Unpacking and Inspection	12
	Matress System Set-up	12
	Replacement Mattress System.....	12
11	Patient Transportation	13
12	Cleaning Procedure	13
13	Care and Storage	14
14	Troubleshooting Guide	15
	Preventive Maintenance	15
15	Disposal	16
	Environment Protection	16
	Disposal	16
16	Warranty	16

1 Symbols

Symbols and Labels on the Product

	Read user manual
	Attention
	CE mark of conformity with EU regulation
	Applied parts type B
	Ingress protection
	Caution risk of electric shock

	<p>The Product identification label is located on the back of the SCU (System Control Unit). The serial number and the model number can be found on the type plate. This information is required when contacting LINET.</p>
--	---

Other Symbols

Structure of instructions:

- ❖ Perform this step.

Result, if necessary.

Structure of bulleted lists:

- List level 1
- List level 2

Warning Notices

Warning words are used in this manual and signify hazards and unsafe practices which could result in personal injury or property damage. See the definitions below for information relating to each word.

⚠ NOTE / IMPORTANT

Indicates a hazardous situation that could result in damage to property if it is not avoided.

⚠ DANGER

Danger indicates an imminently hazardous situation, which, if not avoided, will result in serious injury or death.

⚠ CAUTION

Caution indicates a potentially hazardous situation which, if not avoided, may result in property damage or minor injury or both.

⚠ WARNING

Warning indicates an imminently hazardous situation which, if not avoided, will result in serious injury or death.

⚠ RISK OF ELECTRICAL SHOCK

Do not remove control unit cover.

(AP) EXPLOSION HAZARD

Do not use the control unit in the presence of flammable anesthetics or in the proximity of an oxygen tent.

There is no known risk of adverse effects on the LINET control unit/pump caused by other electromagnetic devices, present at the time of treatment, or vice-versa.

- Refer servicing to qualified service personnel.
- Never drop or insert objects into any opening of the control unit.
- DO NOT SMOKE while using this product and do not use in the presence of smoking materials or open flame. Smoking by visitors in the room will contaminate the system. Therefore, visitor smoking is NOT permitted. Air flowing through air mattress will support combustion. Failure to observe this warning can result in severe fire, property damage and cause injury or death.
- Entrapment may occur. Patient entrapment with bed side rails and mattress may cause injury or death. Mattress MUST fit bed frame and side rails snugly to prevent patient entrapment. Follow the manufacturer's instructions and monitor patient frequently. Proper patient assessment and monitoring, and proper maintenance and use of equipment is required to reduce the risk of entrapment. Variations in bed rail dimensions, and mattress thickness, size, or density could increase the risk of entrapment. Refer to the owner's manual for beds and rails for additional product and safety information. Mattress MUST fit bed frame and bed rails snugly to reduce the risk of entrapment.
- To avoid risk of electric shock, this equipment must only be connected to a supply main with a protecti-

ve earth using the power cord provided with the product.

- Do not heat, steam autoclave, immerse the control unit in liquids or spill liquids or food on or into the control unit. In the event of any spillage, immediately turn off the control unit and disconnect it from the power source. Return the control unit for servicing to a Linet authorised service centre.
- Care should be taken such that the power cord of the control unit is not pinched and does not have any objects placed upon it. Make certain the unit and power cord is not located where it can be stepped on or tripped over.
- Do not modify this equipment without authorisation from the manufacturer.
- Not for use in oxygen-rich environments.
- Before opening or cleaning the control unit enclosure, make certain that the unit is turned off and unplugged from its power source. The control unit enclosure should only be opened by factory authorized qualified technical service personnel.
- Please read this manual before operating any of LINET's Air Therapy systems. If you are unable to understand the manual, please contact your dealer or the manufacturer before attempting to use this equipment. Otherwise personal injury or property damage may result.
- When installing the LINET mattress system, do not exceed the manufacturers rated weight for the mattress or the bed frame. See the bed frame manufacturer's manual for bed frame weight rating.

FCC Statements:

This device complies with FCC part 15.231 for license exempt radio apparatus. Operation is subject to the following two conditions:

- This device may not cause harmful interference, and
- This device must accept harmful interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

2 Contraindications

All air mattress systems are contraindicated for patients with cervical traction or unstable:

- Spinal fractures
- Spinal cord injury

Also:

- Fractures at risk of complication by a moving support surface.
- Trauma patients where spinal injuries have not been excluded or 'cleared'.

CAUTION

For patients undergoing cervical traction the use of the LAL mode is contraindicated and the use of LAL and PT modes should only be under the direction of a suitably qualified physician.

Before placing a patient on an air mattress, a risk assessment by a suitably qualified person should be completed, to ensure the appropriateness of the support provided.

3 Technical Specifications

Electrical Specifications

NOTE

The control unit Power Inlet is used as the means of isolating the equipment from the supply mains on all poles simultaneously.

Parameter	Value
Input Voltage AC	90 ~ 240 VAC
Input Frequency	60 / 50 Hz
Maximum Power Consumption	180 W ± 30 W
Circuit Protection	Dual fused, 250V, 5A Fast blow fuse(s), std. fuses
Fuse Type	Bussmann S500-5-R
Breaking Capacity	@125 VAC is 10kA @250 VAC is 200A
Mode of Operation	Continuous

Performance Specifications

Parameter	Value
Weight capacity	Standard replacement mattress
Pressure zones	455 kg
Max. Flow	1275 LPM / 45 CFM
Support Surface Inflation Time	Within 120 seconds (Max Inflate)

Patient Comfort Pressures / Alternating Pressure

Parameter	Value
Soft Pressure	7 ± 5 mmHg
Firm Pressure	32 ± 5 mmHg

Patient Contact

Control unit and the mattress are **lead free, mercury free** and **latex free** components.

Dartex top sheet is Halogen-free (bromide-free).

Mechanical Specifications

Parameter	Value
Dimensions LxWxH	200 x 87 x 20,5 cm
ProphyCair	200 x 99 x 25 cm
ProphyCair +	200 x 87 x 18 cm
ProphyCair Lo	200 x 99 x 18 cm
ProphyCair Lo +	
Weight	4 kg (CU), 7 kg (mattress)
Power cord	600 cm
Connection	1.27 cm flow single quick disconnect connector
Packaging	1 piece per box

PLEASE CLEAN FILTER EVERY 3 MONTHS OR WHENEVER DIRTY

Remove 2 thumb screws and separate filter foam from filter cover. Wash filter foam using soap and water. Dry and replace filter and filter cover on the unit and fasten screws.

Support Surface

Air cushions are high frequency heat sealed, liquid proof and washable. Base is liquid proof and washable. Top Sheet is low friction, low shear, low force producing, breathable, liquid resistant and highly vapor permeable. Material is 70 Denier nylon or Dartex and a bottom foam or air pad on the bottom of the air cushions enclosed in a cover.

Environmental Specifications

Operating conditions

Parameter	Value
Ambient Temperature	40° to 104° F / 10° to 40° C
Relative Humidity	30% to 75% Non-Condensing
Atmospheric Pressure	700 hPa to 1060 hPa
Protection Against Harmful Ingress Of Liquids	IP3X

Storage and shipping conditions

Parameter	Value
Ambient Temperature	-40° to 158° F / -40° to 70° C
Relative Humidity	10% to 100%
Atmospheric Pressure	500 hPa to 1060 hPa
Protection Against Harmful Ingress Of Liquids	IP3X

Disposal Requirements

This equipment should be disposed of at your local Recycling Centre (Non-hazardous waste) when it has reached the end of its service life. SCU is to be disposed of in accordance with WEEE regulations

FDA Registration

FDA registered company as a manufacturer and as a contract manufacturer. LINET'S quality system meets the requirements of:

- FDA 21 CFR, PART 820- QSR-Current Good Manufacturing Practices (CGMP) for medical devices.
- Products are FDA listed.
- Health Canada Medical Device Licensed.
- ISO 13485 certified company.
- Medicare coded (SADMERC).

Patents

US & International Patents Pending

Certificate of Origin

All products are manufactured in Corona, CA, USA.

4 Safety Instructions

- ❖ Before operating your ProphyCair series control unit, be certain the AC power available at your location matches the power requirements printed on the product identification label on the back of the unit.
- ❖ To avoid electric shock, always plug the power cord of the control unit into a properly grounded power source.
- ❖ Follow the instructions carefully.
- ❖ Ensure any user has read and understood this manual completely before operating the product.
- ❖ Use the mattress system only as specified in this manual and in perfect working order.
- ❖ Use the mattress system only with the correct mains supply (see Electrical Specifications).
- ❖ Use the mattress system only in its original state and do not modify it in any way.
- ❖ Have the mattress system used only by or under supervision of trained and qualified nursing personnel.
- ❖ Replace any damaged parts immediately with original spare parts only.
- ❖ Have the mattress system serviced and installed only by qualified personnel trained and authorized by the manufacturer.
- ❖ Do not exceed the maximum patient weight limit (see Mechanical Specifications).
- ❖ Do not use the SCU near flammable gases. This does not apply to oxygen cylinders.
- ❖ Never cover the SCU while in use.
- ❖ Never handle the mains plug with wet hands.
- ❖ Disconnect the product from the mains only by pulling the mains plug. When pulling the mains plug, always hold the plug, not the cable.
- ❖ Position the mains cable so that there are no loops or kinks in the cable; protect the cable from mechanical wear and tear.
- ❖ Improper handling of mains cable can cause an electric shock hazard, other serious injuries or damage to the mattress replacement system.
- ❖ Do not place the patient directly on the mattress without the top sheet.
- ❖ WARNING: Before opening the control unit enclosure, make sure the control unit is turned off and unplugged from its power source. The control unit enclosure should only be opened by a factory authorized qualified service technician.
- ❖ Smoking by the patient or anyone else around or on system is prohibited. ProphyCair series system uses room air for circulation through the mattress. Smoking will contaminate the system.

5 Standard and Regulations

The mattress replacement system complies with the following standards:

- EC directive 93/42/EEC for medical devices
- RoHS directive 2011/65/EU

The manufacturer adheres to a certified quality management system in compliance with the following standards:

- EN 60601-1:2006/A1:2013
- EN 60601-1-2:2007
- EN 60601-1-6:2010
- EN ISO 10993-5:2009
- EN ISO 9001:2015
- EN ISO 13485:2012
- EN ISO 14001:2015

6 Scope of Delivery

Scope of Delivery

Delivery:

- Upon receipt, check that the shipment is complete as specified on the delivery note.
- Notify the carrier and supplier of any deficiencies or damages immediately as well as in writing or make a note on the delivery note.

List of Parts

- Mattress with Dartex cover - Applied part type B
- SCU (System Control Unit) - Applied part type B
- Mains cable
- User Manual

7 Use

Preparing the Bed for the Patient

WARNING

Risk of injury when putting patient into bed!

Before putting the patient into bed:

- Ensure that mattress is completely and correctly inflated.
- Ensure that mattress is correctly secured with safety straps.

WARNING

Danger of suffocation due to air-impermeable mattress cover!

- Use mattress cover correctly.
- The nursing staff are responsible for the safety of the patient on the mattress cover.

CAUTION

Material damage due to dampness or contamination!

- Ensure that no moisture gets into the mattress.
- Ensure that no food or drink gets into the mattress.
- If moisture, food or drink get into the mattress, notify LINET® Service.

Putting the Patient into the Bed

Skin Assessment

Before placing a patient on any air mattress, a skin assessment by a qualified person should be completed to ensure the appropriate support surface is used. Always follow local facility protocol for skin and risk assessments.

- Lay patient on mattress.

For an ideal lying position:


- If additional blankets or sheets are used, make sure that ease of movement is sufficient.
- Ensure that blankets, sheets, clothing etc. do not cause pressure sores.

⚠ WARNING


Due to the many types of bed frames and side rails available, consideration must be given to the safe and appropriate use of side rails and mattress heights when using mattress replacements and overlay systems. A risk assessment should be carried out in accordance with your local instructions.

8 Product Description - Control Unit

These Air Therapy Systems are designed to provide continuous air therapy support at required patient comfort levels and are used to inflate an overlay mattress or replacement mattress system. The Thermoplastic 94V-0 fire retardant enclosure houses a blower, main PCB, display panel, a short circuit /over voltage protection with dual fuses, a quick connection and disconnection (CPR deflation) coupling connector (**F**) with 1.27 cm flow port with RTPM (real time pressure monitoring) technology sensors, and a 19.7ft (6 m) long detachable hospital grade power cord.

Controls and Indicators

Operating Instructions

POWER (number 1)

- ❖ During initial power up when power cord (1) is plugged into the power source, the control unit will display illuminated amber light indicating that the unit is in Standby mode.
- ❖ Press button ON/OFF Power key, the (amber - OFF) Standby LED will turn off and the (green - ON) LED turns on. The blower will turn on.

PATIENT COMFORT CONTROL LEVEL (number 2)

- ❖ The system is designed for patients weighing between 15 Kg. ~ 455 kg. The pressure level setting is governed by the patients weight. The scale on the right hand side of display label suggests the level setting appropriate to different weights, and should be referred to in the first instance when setting up the mattress for use.
- ❖ Pressing the Comfort Control Keys (2) allows the patient comfort level to be adjusted."-“ reduces the pressure setting level, “+” increases the pressure setting level.
- ❖ Depending on the desired patient comfort level, the controller will increase or decrease the speed of the air blower to provide the appropriate air flow into the mattress to maintain the desired pressure in the mattress.
- ❖ Once the mattress has been fully inflated using MAX inflate mode, lay the patient on the mattress and set the Comfort Control Key(s)to the desired patient comfort level. Wait 2 minutes for the mattress pressure to stabilize and verify the appropriate pressure required to support the patient by performing a simple “four finger check”. Make sure that the patient is lying flat on their back in the centre of the mattress. Place four fingers between the air cushions directly underneath the sacral region of the patients body. There should be a minimum of 3 to 4 finger width clearance between the bottom of the patient and the safety foam base.
- ❖ Repeat this procedure until the desired patient comfort pressure is achieved.

DISPLAY (number 3)

- ❖ Indication of the pressure setting selected using +/- buttons (2)

POWER FAIL (number 4)

- ❖ In the event of a power loss, the microprocessor will activate an audio-visual signal to alert the caregiver by flashing the amber “POWER FAIL” LED and turning on the buzzer. Once the power is restored to the control unit the audio-visual signal will cease and unit resumes operating its set mode. When deliberately disconnecting the control unit the Alarm Silence (8) should be used to cancel the audio alarm.

FOWLER (number 5)

- ❖ This mode is selected by pressing the (number 5) Fowler key if the patient is in upright mode / fowler position (when the bed frame is articulated to fowler position of 35° or greater). In upright mode the unit will maintain 80% more pressure than the set pressure (up to max 32 ± 5 mmHg) in the torso section of the mattress in order to keep the patient supported and help prevent the patient bottoming out.

- ❖ AUTO FOWLER. When the head section of the bed frame is articulated to greater than 35°, the unit will automatically enter into fowler boost mode, without the caregiver's or the patient's intervention. Its activation is signified by the illumination of the green LED adjacent to key number 5. In this mode the unit will maintain 80% more pressure than the set pressure (up to max to 32 ± 5 mmHg) in the torso section of the mattress in order to keep the patient supported without the patient being bottomed out. This is achieved by a low power wireless sensor in the mattress. The Auto Fowler is factory programmed to a 35° angle. It uses a +3.6V Lithium battery that can be replaced using LINET part number 100622-S.

Should the Auto Fowler Boost function not be desired, it can be toggled ON/OFF by pressing the Fowler Key (number 5). Also, in the unlikely event that the sensor fails, this mode can be manually selected by pressing the (number 5) Fowler key for when the patient is in any position.

MAX INFLATE (number 6)

- ❖ Press MAX Inflate (6) key, the (green) LED will turn on. This mode is used to rapidly inflate the mattress. During this mode a series of beeps will sound every 3 minutes as a reminder that MAX Inflate mode has been activated. During this mode, the entire mattress will be pressurized to a pressure greater than Level 9. The mattress will inflate to its normal size within 120 seconds. MAX Inflate mode will deactivate after 15 minutes.
 - For rapid inflation of the mattress, press (number 6) "MAX Inflate" key (green LED turns on).
 - For extra firm support during patient ingress or egress, patient wound care, patient turning or patient cleaning it is recommended to set the mattress pressure to maximum pressure by pressing (number 6) "MAX Inflate" key.

PROGRAM (number 7)

- ❖ LOW AIR LOSS (LAL)
 - To set LAL mode, press Select key (number 7) until the "THERAPY" LED (LAL) lights up.
 - In LAL mode all the air cushions in the mattress will be maintained at the same pressure and can be adjusted to the desired patient comfort control level.
- ❖ PULSE (PT)
 - To set PULSE mode, press the Select key (number 7) until the "PULSE" LED (PT) lights up.
 - During Pulse mode, the pressure in each air cell will be decreased by 20% for 5 seconds and then increase to the prior pressure setting for 30 seconds.

ALARM SILENCE (number 8)

- ❖ Mutes the audio alarm. The green indicator will come on to remind the user that an audio alarm has been muted. This indicator will go off if the mute is canceled by pressing the key a second time.

LOCK OUT (number 9)

- ❖ Control unit functions (including power switch) can be completely locked out to prevent tampering of settings by simply pressing and holding the Lock key until the light comes on (approximately 5 seconds). To deactivate, simply press and hold the lock key until the light turns off. Activation and deactivation will be confirmed by a double audio tone.
- ❖ Max Inflate mode is still functional in Lock Out mode.

LOW PRESSURE (number 10)

In the event of hose disconnection or other problem that causes major air leakage, this will activate an audio-visual signal to alert the caregiver by flashing the amber "LOW PRESSURE" LED and turning on the audio alarm. Once the low-pressure problem is fixed, the audio-visual signal will cease and the unit will resume its operating set mode. The audio alarm can be muted using the Alarm silence key.


CPR

- ❖ To deflate the mattress for a CPR procedure, hold down the grey connector latch and pull the hose from the control unit connector. The audio and low pressure visual alarms will start. If required the audio alarm can be muted using the Alarm silence key.

9 Product Description - Mattress

Self-contained mattress replacement system/ mattress overlay system (A) with easily detachable components for cleaning. Flame retardant, water repellent, mildew resistant, low friction and low shear, 5" or 8" high (inflated) detachable lateral tubular air cushions (B) (16~22). Detachable zippered or strapped highly breathable, reusable, quilted Top Sheet composed of two layers; Top layer is made out of urethane coated 70 Denier nylon or Dartex material. Bottom layer is made out of quilted material. It is flame retardant, water repellent, vapor permeable, anti-microbial, low friction and low shear (C). 2" safety air pad or convoluted foam pad enclosed in a foam/air pad base (D) to support the patient in the event of loss of air pressure in the mattress. The mattress has a hose assembly (E) with easy to use quick connect and disconnect connector coupling (F).

Mattress Weight Capacity

NOTE: See bed frame manufacturer's manual for bed frame weight rating.

ProphyCair systems are used for patients who weigh between 15 kg and 455 kg.

WARNING

It is highly recommended that the ProphyCair system always be installed on medical bed frames that are equipped with standard hospital side rails or assist rails. Please raise all side rails on the bed and lock them in position after the patient is on the mattress. Health care professionals assigned to each case should make the final determination whether side or assist rails are warranted after assessing patient risks of entrapment and falls in accordance with State patient restraint legislation or facility interpretation of such legislation.

Check that all air hoses and power cord are clear of moving bed components before placing a patient on the bed and that the mattress system is fully inflated. Operate all bed frame motorized functions through their full range of motion to be certain that there is no pulling, interference or pinching.

Mattress MUST fit bed frame and side rails snugly to prevent patient entrapment.

10 Unpacking the System

NOTE

When opening the large system box or the small control unit box, ensure that the object used to open the box does not penetrate and damage the components inside.

Components Supplied

- Complete Replacement Mattress System Box
- 1 Control Unit box & 1 Mattress

Unpacking and Inspection

- ❖ Before accepting and signing for your shipment, please inspect the box or boxes for external and internal damage. Verify that the number of boxes listed on the packing list matches the number of boxes received. Verify that no components in your shipment are damaged or missing. Report any missing boxes, components and or damage to the transportation carrier immediately.

Matress System Set-up

NOTE

Make sure that the hose end of the mattress is towards the foot of the bed.

NOTE

When installing the LINET mattress system, do not exceed the manufacturers rated weight of the mattress or the bed frame. See the bed manufacturer's manual for bed frame weight rating.

NOTE

Care should be taken such that the air inlet vent on the control unit is not covered, and the control unit is not placed on the floor in such a manner that it is a hazard to caregivers, patients or the public or lowering of bed frame.

NOTE

Before using the ProphyCair Mattress System, please remove current mattress, mattress replacement systems or overlay systems from the bed.

CAUTION

Ensure that the mattress is fully inflated in MAX inflate mode before placing the patient into bed.

Replacement Mattress System

- ❖ When using the replacement mattress system, the mattress is placed directly on the bed frame and remove any existing mattress.
- ❖ Unroll the Replacement Mattress and place it on the bed frame (**BF**).
- ❖ There are ten fixing black straps with buckles (**G**), two straps at the head of the mattress, two on the foot of the mattress, and three on each side of the mattress. Loop each strap around the bed deck and secure using the buckle (trident clip). Note: Make sure Head, Knee and Foot sections of the bed can be raised and the straps are secured to the bed deck and not to the frame. Once the mattress is strapped, tuck the exposed straps under the mattress.
- ❖ Pull out the hanging hooks on the back of the control unit (**H**) and suspend the control unit from the footboard (**FB**) of the bed (**BF**). If the bed you are using does not have a footboard, place the control unit (**H**) on its base (not on its back where the filter is located) on a flat surface in front of the bed near the foot of the bed frame (**BF**).
- ❖ Uncoil the power cord (**I**) and plug the cord into the appropriate AC power source, which is properly grounded. Plug the other end of the power cord into the control unit and press it in place.

NOTE

Care should be taken such that the power cord of the control unit is not pinched, or has any objects placed on it. Ensure it is not located where it can be stepped on or tripped over. Make sure the control units power inlet connection is positioned to easily disconnect the power cord from the unit.

- ❖ Connect the mating connector (**F**) on the mattress hose assembly (**E**) into the insert on the control unit connector and lock it in place.

11 Patient Transportation

Air Mattress

- To transport a patient without removing the patient off the bed, turn off the control unit, disconnect the power cord from the power source and roll it up on the control unit securely.
- The mattress will deflate within a few seconds after the control unit has been turned off.
- The mattress has a 2" foam or air pad to provide support to the patient when the mattress is deflated. It is not recommended to keep the patient on the mattress for long periods of time when the mattress air cushions are deflated.

12 Cleaning Procedure

Mattress cleaning

All components such as manifold, hose assembly, air cushions, top sheet, and foam base are interchangeable and can be easily cleaned or detached for laundry.

WARNING:

Before attempting to clean the control unit, turn off and disconnect the power cord from the power source.

WARNING:

DO NOT HEAT, STEAM AUTOCLAVE, OR IMMERSE THE CONTROL UNIT IN LIQUIDS

Control Unit:

- ❖ Wear eye goggles and protective gloves before starting the cleaning procedure.
- ❖ The following germicidal detergents / disinfectants are recommended by the EPA as hospital disinfectants.

Johnson Wax, Virex 128, EPA Registration Number 47371-130-4822.

*Quaternary Detergent-Disinfectant by
Airkem Professional Products, Division of
Ecolab, Inc., Ecolab Center, St. Paul,
Minnesota. EPA registration number: EPA
42964-5.*

*Hi-Tor Germicidal Detergent by Huntington
Laboratories, Inc. Huntington, Indiana.
EPA registration number: EPA # 303-91.*

NOTE:

A spray bottle of fresh disinfectant / detergent solution should be prepared daily to clean the control unit.

- ❖ Prepare the required amount of solution by following the preparation instructions provided with the germicidal detergent /disinfectant solution.
- ❖ Pour required amount into a spray bottle.
- ❖ Use a brush or cloth to wipe off dust. If necessary, spray the exterior of the top and the bottom enclosures, power cord and the cord plug with the prepared disinfectant / detergent solution. Using a damp cloth, wipe down the sprayed surface cleanly. Note: Do not spray excess amount of solution on the control unit.
- ❖ Once the control unit is clean, wipe the unit, power cord, cord receptacle, and the cord plug with a clean dry cloth.
- ❖ Place the control unit to dry in a cool, dry area for an hour before operating or storing the unit. If the control unit is not used immediately, place the control unit in a plastic bag and store it in a storage area.
- ❖ After the cleaning operations are completed, remove and dispose the protective gloves appropriately. Wash your hands thoroughly with antibacterial soap.

Air Mattress:**⚠ WARNING:**

To prevent damage, please remove the Auto Fowler sensor from the base of the mattress before laundering. The Auto Fowler sensor is located in a small pocket inside the base of the mattress located on the patient's right hand side, 30 - 38 cm from the head end of the mattress.

- ❖ Wear eye goggles and protective gloves before starting the cleaning procedure.
- ❖ Follow steps 2 through 4 above to prepare disinfectant solution.
- ❖ Using a damp cloth, wipe down the air cushions and the mattress base. Once the air cushions and the base are clean, wipe them down with a dry cloth.
- ❖ Air cushions should be washed periodically; top sheet will require more frequent washing. Set wash cycle to "heavy load" with warm water. Once the water is full, add manufacturer- suggested quantity of laundry detergent and/ or standard hospital disinfectants. If the air cushions or the top sheet becomes soiled with human waste or blood, clean immediately by wiping down. Use hospital recommended laundry detergent and/ or disinfectant per manufacturer's instructions. Note: Use non-chlorine bleach detergent.
- ❖ Shake cushions gently to remove excess water from inside the air cushions. Dry the cushions/top sheet on the lowest settings on the dryer until completely dry.
- ❖ Leave the mattress to dry in a cool, dry area for an hour before using or storing. If the mattress is not used immediately, roll the mattress and insert it into a plastic bag and store it in a storage area.
- ❖ After the cleaning operations are completed, remove and dispose the protective gloves appropriately. Wash your hands thoroughly with antibacterial soap.

13 Care and Storage

- ❖ When the control unit is not in use, turn off the unit, disconnect the power cord from the power source and wrap the cord around the control unit. Place the control unit and the power cord in a plastic bag and cable tie it to keep the unit dust-free.
- ❖ Fold or roll the previously dried air mattress and place the mattress in a plastic bag. Cable tie the plastic bag to keep the mattress dust free.
- ❖ Store the control unit in a storage area designated for medical electronic product storage.

14 Troubleshooting Guide

⚠ THE FOLLOWING INFORMATION IS FOR AUTHORIZED SERVICE FACILITIES AND QUALIFIED SERVICE PERSONNEL ONLY.

LINET can provide technical support to Linet qualified technical personnel. Contact LINET service department for more information.

PROBLEM	CAUSE	SOLUTION
A. Mattress not inflating properly	1. Mattress hose disconnected 2. Air hose kinked or split 3. Major leak in the air cushion or overlay pad 4. Kinked or split manifold 5. Control unit not working 6. Blower malfunction	1. Connect hose connectors and lock them in place 2. Unkink hose or <i>replace split hose</i> 3. <i>Replace leaking air cushion or overlay pad</i> 4. Unkink manifold or <i>replace split manifold</i> 5. Send control unit back to factory for repair 6. Send control unit back to factory for repair
B. No Power	1. Control Unit OFF 2. Power cord disconnected 3. No power in the power source 4. Power outage 5. Blown fuse	1. Check power source and turn unit on 2. Connect cord to power source 3. Check power source has power and turn it "ON" 4. Wait till the power source has power 5. Send control unit back to factory for repair
C. Control Unit not responding	1. Unit locks up	1. Unplug power cord from control unit, wait 15 minutes and plug cord back to reset the control unit

Preventive Maintenance

It is important to periodically test the control unit to verify its functionality. If the units air pressure reading is out of specification it can result in poor or reduced patient support. Contact Linet Service for training on how to perform these tests.

⚠ PLEASE CLEAN FILTER EVERY 3 MONTHS OR WHENEVER DIRTY

Remove the 2 filter screws from the back of the unit and separate filter foam. Wash filter foam using soap and water. Dry and replace filter back on the unit and fasten screws.

⚠ NOTE:

All preventive maintenance service, performance and electrical tests, or repairs should be performed only by authorized and qualified technical personnel.

15 Disposal

Environment Protection

LINET® is aware of the important role that the protection of our environment plays for future generations. The materials of this product are environmentally compatible. It does not contain hazardous substances on the basis of cadmium, mercury, asbestos, PCB or CFC. The noise emission and the vibrations meet the directives for healthcare premises. The packaging materials are produced according to the respective directives. Dispose of the packaging material according to the symbols and by delivering it to an authorised person.

Disposal

To dispose of the appliance (SCU): When you dispose of your appliance (SCU) do not put it into the household waste.

Send the appliance (SCU) to the recycling of electrical appliances.

16 Warranty

LINET® will only guarantee the quality and reliability of products, that are periodically maintained and used in accordance with safety instruction.

If there occurs the serious defect, which is not possible to be repaired within the maintenance:
Do not continue to use the mattress system.

This product is covered by a 24-month warranty. Beginning of warranty time is understood the date of product's dispatch from company LINET® to the final user. The warranty covers all material and manufacturing-related failures and errors. Failures and errors caused by incorrect use and external effects are not covered. Justified complaints will be fixed free of charge during the warranty period. Proof of purchase, with the date of purchase, is required for all warranty service. Our standard terms and conditions apply.

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