



sella

Multifunctional Chair





Multifunctionality, mobility and mobilization


sella - the solution for diverse areas in clinic and nursing care

Patients and residents are often very limited in their mobility.

The reasons for this situation are manifold, such as prolonged illness – combined with a stay in hospital or advanced age. In many cases, the consequences are reduced strength and limited mobility, which can often lead to an insecure gait with an associated risk of falling.

With early mobilization attempts and movement-enhancing measures in the clinic and nursing care, independence and self-activity shall be promoted, but also consequential damage – such as increasing immobility with the risk of contractures and decubital ulcers – shall be significantly reduced. A first step towards promoting mobility may entail a transfer from a hospital/care bed to an adequate multifunctional chair.

The sella multifunctional chair is equipped with functions that promote mobility and mobilization and provide optimum support for both patients/home residents and caregivers in terms of safety, ergonomics and comfort. In addition, sella's excellent driving characteristics make it ideally suited for a wide range of transport/transfer situations within the hospital and/or care facility.




wissner-bosserhoff | sella

Flexibility as needed and application



- Comfortable, safe entry and exit height fall prevention
- Passive and active safety elements: side guards, armrests, pillows, storage table, central brake and foot pedal
- Positioning as required: sitting, lying and resting



Flexible table - strengthening the seating position and promoting independence and self-reliance



Side protection, padded safe sitting and lying comfort



Workflow Optimization

- Mobility and mobilization in one product
- Extensive accessories for clinic and care
- Support of mobilization steps under ergonomic aspects lifting hydraulics for height and inclination of the seat



Comfortable resting and sleeping positioning or support with transferor examination situations



Lifting hydraulics for height adjustment with seat inclination for a safe entry and exit of the patient as well as back-friendly care support.



Process Optimization



- Flexible areas of application in clinics and nursing homes
 - Extensive model variants for indoor and outdoor use
 - Easy and safe maneuvering in different environments



Outdoor: Maximum comfort for long distances or for uneven surfaces thanks to large casters



Indoor: Safe and easy maneuvering in the smallest space thanks to small castors at the front

Hygiene



- Removable, antimicrobial upholstery
- Infection control through easy-to-clean surfaces
 - · Easy removal of the seat upholstery



Infection control through easy cleaning of the upholstery elements

wissner-bosserhoff | sella wissner-bosserhoff | sella

Portfolio

Operation area

Drivability

Modell

Maximum patient weight

Dimensions

Outdoor Indoor SLA-AX SLA-DX SLA-DX Lift carriage with rotating front castors Fixed carriage with big front castors 2 big front castors -4 small castors for outdoor transfer for the optimal use in small rooms 150 kg 150 kg Length in basic position 910 mm 1610 mm Length in lying position Width 685/705 mm* Maximum lift (of height adjustment) 205 mm 490-505 mm** Minimum height 490-745 mm** Maximal height Back rest adjustment angle 12° bis 70°

Colour options standard

(Accurate colouring is only possible by means of decorative samples.)



Baltic (008278) Synthetic leather



Apricot (008294) Synthetic leather



Footrest adjustment angle

Arm rest adjustment angle

Creme (008291) Synthetic leather



Pink (8048) Fabric



12° bis 87°

-25°, 0°, 25°, 90°

Perweiß (6017) Fabric

Special colours

Other colours are available for an additional charge. (Exact colouring is only possible with decor samples). The following special colours are shown in the brochure:



Blue



·

Accessories



Table



Table with safety barrier



Table safety lock



Table holder



Telescopic infusion pole



Holder for infusion pole



Stacking pocket



Urine bag holder



Headrest, big



Side support (pair)




Independent adjustment of backrest



Independent adjustment of legrest





Member of LINET Group

wissner-bosserhoff GmbH

Hauptstr. 4-6 | D-58739 Wickede (Ruhr) | Germany

Phone: +49 2377.784-0 | Fax: -163 | info@wi-bo.de | wi-bo.com

