

Special beds

tom

Children's bed

- Mattress dimensions 132 x 67.5 cm and backrest adjustment up to about 60°
- Height-adjustable and tiltable mesh-base patient surface (Anti-Trendelenburg / Trendelenburg position)
- Initial fault tolerant side guards, adjustable on 3 levels
- Extensive choice of colors and accessories



rooming in

Infant bed

- Transparent and removable infant basket with ventilation holes
- Continuously height-adjustable and tiltable bed surface
- Suitable for moving under a universal hospital bed



sprint

Stretcher

- Excellent movability thanks to integral castors and optional 5th castor
- Two-part, hydraulically adjustable and tiltable patient surface with backrest adjustment
- Optional X-ray-capable patient surface for C-arm or standard X-ray exposure
- One-of-a-kind mattress system with three different variations (standard, comfort, transfer)
- Closed lifting column construction ensures efficient cleaning
- Several wall bumpers and options for attaching accessories
- Safe working load of 230 kg



wissner-
bosserhoff

Functional elegance in nursing care

wissner-bosserhoff GmbH
Hauptstraße 4-6
58739 Wickede (Ruhr), Germany
Phone: +49(0)2377/784-0
Fax: +49(0)2377/784-163
E-Mail: info@wi-bo.de

www.wi-bo.com

wissner-
bosserhoff

Functional elegance in nursing care



Bed portfolio
Hospital

wissner-bosserhoff | Art.-No.: 97-0900 | Edition: 09/2011, Rev.0 | Print: msp druck und medien | Slight deviations in colour are possible. Technical changes reserved.

Universal beds – hotel design

image

Innovative mobilization in rooms with additional services

- Continuous side guards can be completely lowered beneath the patient surface when not in use, for unrestricted access to the patient and accessories and for smooth patient transfer
- Mobi-Stick and Mobi-Lift® (optional) and pre-programmed bed egress position support self-mobilization and physical therapy, for instance in surgical and orthopedic wards
- Prevention of bedsores thanks to double retraction of the backrest and leg rest
- Fall prevention thanks to four side guard positions and a brake alarm (optional)



eleganza smart Hotel

Safety and flexibility

- More than just compliance with the new standard IEC 60601-2-52 thanks to central gap solution smaller than 60 mm, resulting in a high degree of patient safety in all working positions, even without restricting central gap protector
- Divided side guards and integrated operating elements support mobilization
- High degree of modularity
- High abrasion resistance with low-weight side guards



altura

Economical and cozy

- Economical start to a modern and cozy hotel design
- Continuous side guards = 100% protection in every working position
- Reduction of FDM* possible by lowering of side guards at the foot end
- High functionality with useful additional features: 5th castor for easy transport, foot switch for bed operation and removable ABS patient surface for ease of nursing care



latera thema

Lateralization for ease of nursing care

- Fewer physically strenuous activities thanks to lateralization up to 15°
- Time-savings of up to 15% with routine nursing tasks
- Significant reduction in the strain on nursing staff, especially during night and weekend shifts
- High return on investment thanks to sustainable process cost savings and high level of motivation on the part of the nursing staff
- Continuous side guards for seamless safety
- Cozy design with numerous décor variations



Lateral tilting makes daily nursing care easier

It is no secret: Statistically speaking, nursing care is one of the most high-risk occupations. In addition to infections and stress, nursing staff particularly often suffer from back pain and blockages in and deformations of the spinal column, often with permanent damage. >>

Intensive care beds


multicare

The new first-class intensive care bed

- In-bed therapy**
1. Kinetic therapy concept by means of lateral therapy up to 30°
 2. Bedsore therapy system virtuoso
 3. Upright sitting position in the bed thanks to the orthopneic position
- In-bed diagnosis**
4. C-arm compatible backrest
 5. X-ray permeable backrest with integrated X-ray cassette holder
 6. Integrated weighing system with storage of weight data and analysis tool as well as fall prevention thanks to bed egress alarm

- Flexibility**
7. Easier transport due to central 5th castor
 8. Fall prevention thanks to automatic bed brake (i-brake)

- Further advantages of multicare standard**
- Product designed in accordance with the new standard IEC 60601-2-52; in particular in terms of the initial fault tolerance, side guard concept and electrical safety
 - LCD touch screen for easy single-hand operation
 - Safe working load of 250 kg
 - Lateral tilting up to 15° for ease of nursing care



eleganza 3XC

Perfect for intermediate care and intensive care units (ICU)

- Retractable cassette holder on the side, for lung X-rays
- C-arm compatible from the upper body to the pelvis area
- Ergoframe patient surface system for bedsores prevention
- Compatible with active mattress systems
- Mobilization aid Mobi-Lift®
- Integrated weighing system (optional)

Universal beds – clinic design

latera acute

Lateralization for ease of nursing care

- Fewer physically strenuous activities thanks to lateralization up to 15°
- Time-savings of up to 15% with routine nursing tasks
- Significant reduction in the strain on nursing staff, especially during night and weekend shifts
- High return on investment thanks to sustainable process cost savings and high level of motivation on the part of the nursing staff
- ¾ side guards allow for the use of mattress systems with a thickness of up to 17 cm
- Modern design made of high-quality plastic materials



Your 10 advantages of lateralization:

Making nursing easier

1. Changing sheets and dressings
2. Basic care
3. Mobilization
4. Transfer (to the bed, stretcher, transport chair)

Prevention and care

5. Wound care
6. Prevention of bedsores
7. Breathing more easily
8. Burn injuries
9. Care of infants
10. Receiving visitors

>> Back ailments, in particular, can be significantly reduced by lateral tilting, since mobilization tasks are controlled by the nursing staff without physical effort, thereby easing the strain on their backs. The lateral tilting on both sides of the latera and multicare beds makes it easier to carry out nursing tasks recurring on a daily basis. It not only reduces the effort required of nursing staff, but also makes it possible to care for the patient in a respectful and careful manner.

eleganza 3

Optimally suited for all wards, from standard to intermediate care

- Flexible side guard concept allows for
 - Protection adjustment in line with individual risk of falling
 - Mobilization and standing-up aid, in accordance with the recommendation of the expert standard on fall prophylaxis**
 - Complete protection for patient safety and secure attachment in accordance with the BfArM (German Federal Institute for Drugs and Medical Devices) recommendations***
 - ¾ side guard protection does not represent a FDM*
 - Unrestricted and safe access to the patient in the nursing position
- Prevention of bedsores in the pelvis and shoulder area, as well as permanent support of the head in all bed positions thanks to the Ergo-frame® patient surface concept
- Stationary head section and central 5th castor to avoid damaging walls and door frames
- Mobi-Lift® (optional) for autonomous mobilization



eleganza smart Klinik

Safety and flexibility

- More than just compliance with the new standard IEC 60601-2-52 thanks to central gap solution smaller than 60 mm, resulting in a high degree of patient safety in all working positions, even without restricting central gap protector
- Divided side guards and integrated operating elements support mobilization
- High degree of modularity
- High chemical resistance with low-weight side guards



eleganza smart junior

Particularly for young patients

- Kids and youth bed with surface dimensions of 170 x 80 cm
- Integrated ¾ side guards for fall prevention
- Initial fault tolerance thanks to the GO safety button



* FDM – freedom-depriving measures | ** c.f. expert standard fall prophylaxis, dated February 2005, Chapter 2.6.2.7., p. 71/72. Information available at www.dnqp.de
*** BfArM – Information for specialists from the German Federal Institute for Drugs and Medical Devices dated 12/06/2003