

  **wissner-**
  **bossert**

Functional elegance in nursing care

SafeLift

Intuitive simplicity





Photo: René Schwerdtel, model: Friedhelm Adorf *1943 | cover photo (portrait): Jim Rakete, model: Ernst Strätmann *1933

Improving and maintaining mobility

SafeLift | Revolutionary bed operation for residents

A key goal in in-patient care for the elderly is to maintain or improve mobility. More mobile residents are less prone to falling, further secondary diseases such as pressure ulcers are prevented, and their health is kept in optimum condition. This increases their self-esteem and an independent lifestyle is less restricted.

SafeLift offers the simplest operation so that residents can adjust their bed independently and set it to the best mobilization position for them. At the same time, they can support themselves when standing up on the SafeLift controller. This promotes residents' self-mobilization, reduces staff calls for bed adjustment and frees up care.



Simple bed operation for residents



Standing-up and mobilization aid



Cable- and battery-free controller



Customized for each resident

Simple bed operation & safe mobilization

SafeLift | Product benefits at a glance



01

Adjustment buttons for bed height and backrest
only the features which residents really needs

02

Large buttons
easy to feel and fluorescent lighting

03

Firm position on a SafeFree® side guard
make it easier to find

04

GO button protects against
unintentional bed adjustment

05

360° rotation
and flexible adjustment
on the SafeFree® side guard

06

Ergonomic design
suitable for any size of hand and easy
accessibility of the operating buttons

07

Strong stability
provides firm grip when standing up



Mobilization
and standing-up aid

Residents can set their bed
independently to the best position
to get out and support themselves
on the SafeLift controller.



No cable
or battery

Thanks to established radio technology,
the SafeLift controller can be operated
without any fault-prone cables.



Automatic
night light

and fluorescent buttons make it easier
to get around in the dark.




SafeLift – intuitive simplicity

Added value for residents, care and care-home management

 Resident



 Nursing staff



 Management



Mobilization and independence

- Accessible to residents when lying down, individually positionable and usable as a support
- Fall prevention – no tripping hazards from hanging cables
- Position of the SafeLift controller individually adjustable for each resident
- Suitable for right- and left-handed people, as well as people who get out of the right or left side of the bed
- Automatic night light makes it easier to get around in the dark

Promotes relief & mobility in line with expert standards

- Promotes the mobility and independence of the resident
- For routine tasks, always ready to hand for care, too, care positions as usual by hand switch
- Automatic night light
- Simple and tool-free installation on every SafeFree® side guard*
- Extremely easy to adjust to individual resident needs
- Relief on care staff thanks to reduced staff calls regarding bed adjustments

Flexibility and economic efficiency

- Low follow-up costs – no fault-prone cables and no battery
- SafeLift controller can be attached without tools to any SafeFree® side guard
- Can be used on any bed
- Compatible with the SafeSense® Bed Exit Assistance System
- Can be retrofitted to all wissner-bosserhoff care beds with SafeFree® side guards*



Members of LINET Group

wissner-bosserhoff GmbH

Hauptstraße 4-6 | 58739 Wickede (Ruhr) | Germany

Phone: +49 2377.784-0 | Fax: -163 | info@wi-bo.de | wi-bo.com

