

■ ■ wissner-
■ ■ bossenhoff

Functional elegance in nursing care



Accessories

Nursing care





Accessories

With love for detail

When it comes to care for the elderly, a good human interaction between nursing staff and residents is very important. wissner-bosserhoff understood that an open-minded and friendly personal interaction can be supported by helpful and innovative products. When less time is needed for necessary jobs and tasks, more time is left to care for residents. We have developed useful accessories for our products, which contribute to better quality and efficiency in the field of nursing care.

Contents

Lights	4 - 5
Bed exit system	6 - 7
Mobilization	8 - 11
Storage	12 - 13
Protection	14 - 17
Care	18 - 19

Lights

The undulating shape of the lampshade, with its partially translucent character, is prominent as a special design element of soluna lights. Equipped with a swivel arm, the light allows the light focus to be precisely positioned close to the object to be illuminated, regardless of whether the light is attached to the bed or the bedside cabinet. Either halogen or LED illuminants can be chosen as the light source. Advantage of LEDs: a long life, low energy consumption, no heat generated and absolutely safe operation.



50-0789
Design reading light soluna
HALOGEN, incl. 12 V transformer

50-0773
Design reading light soluna
LED, incl. power supply unit

50-0550
Design reading light soluna
LED, without power supply unit



50-0677
LED reading light lymera
(with and without power supply unit)

Transparent lampshade
(Variant -0000)
Brighter lampshade
(Variant -0001)



01-010812, 01-010813
Switchable underbed light
(Night light function)



Safety by day and night

SafeSense® – the intelligent bed exit sensor system



wissner-bosserhoff, your reliable partner for the reduction of freedom-depriving measures (FDM): With SafeSense®, wissner-bosserhoff presents an intelligent assistance sensor system, which reliably helps reducing freedom-depriving measures day and night. It is equipped with an automatic night light, which is activated when the resident leaves the bed and is deactivated when the resident is back in bed again. SafeSense® is compatible with standard nurse call systems and can be retrofitted to wissner-bosserhoff nursing home beds manufactured in 2002 or later.



Automatic night light

- Night light is automatically switched on when the resident leaves the bed
- Night light is switched off again as soon as the resident is back in bed

Individual alarm system

- Freely adjustable immediate alarm from 0 seconds to 30 minutes
- Adaptation to the resident's mobility

Integrated or available as retrofit kit*

- Accessory option for current WIBO nursing home beds
- Retrofit kit for electrical WIBO nursing home beds manufactured in 2002 or later
- Available with cable or wireless




**More safety –
less freedom-depriving
measures (FDM)**

- Faster response time in the event of falls
- Greater safety during night and weekend shifts
- Fewer inspections
- Ideal in combination with sentida low nursing beds and SafeFree® side guards

SafeSense® bed exit system

SafeSense® Sensor mat 2.1




02-000876
SafeSense® sensor mat 2.1
80 x 20 cm



04-1851
Digital bed exit box 2
cable or wireless version

Also available as a set:
Sensor Mat 2.1 & Bed Exit Box 2
10-0871
Cable version
10-0872
Radio version

SafeSense® Connection kit




04-1508
Adapter cable with on/off switch and suitable plug for connection to the nurse call system (customer-specific) for the cable and wireless version




04-1608
Adapter cable without on/off switch with suitable plug to connect to the nurse call system (customer-specific) for the cable and wireless version

SafeSense® Retaining plate



01-011193
Retaining plate for DU09 motor



01-011194
Retaining plate for Q7 motor

Mobilizing

Independent mobilization reduces the risk of falling. Our accessories help the resident to maintain their independence. We would all like to remain independently mobile for as long as possible. This is why it is so important to preserve and promote the individual mobility of care-dependent residents in line with the standard, developed by experts, known as „Preservation and Promotion of Mobility in the Field of Care“. The benefits for resident range from social participation to health promotion.




10-0799
SafeLift Mobilisations-Assistent,
Set including box

SafeLift bed operation




10-0811
SafeLift Controller




10-0808
SafeLift Box

SafeLift Retaining plate



01-011193
Retaining plate
for DUO9 motor



01-011194
Retaining plate
for Q7 motor

Revolutionary bed operation for the resident

SafeLift | Small handle with a difference




- Maintain & improve mobility
- Simple bed operation
- No cable, no battery
- An automatic night light
- Easy installation, retrofitting

Mobilizing



10-0586
Serving tray with bookholding rails



10-0715
Standing-up aid for 09 and 10 side guards, as well as for ¾ side guards





10-0728
Accessory adapter for standing-up aid (10-0715)



04-0674 / 10-0502
Lifting pole with oval trapeze handle and belt retractor



Quick and convenient access for personnel and residents



10-0422 (A)
Accessory holder to attach to continuous (09) side guards
10-0421 (B)
Handset holder can be plugged into the side guards accessory holder (10-0422)



10-0419 (A)
Adapter to attach to the divided (10) side guards and handset holder SD
10-0421 (B)
Handset holder can be plugged into the side guards accessory holder (10-0419)

Storage

An extensive range of accessories to supplement our products and to support everyday nursing care completes the overall portfolio. This makes nursing care considerably easier for the staff and valuable time is saved, which can then be used to increase the quality of care.



10-0617
Integrated linen holder,
compatible with all beds



Integrated linen holder, pulled out



10-0719
Linen holder compatible with all bed types



Protector to close the foot-end gap of 220 cm beds
10-0580 for extension via cable pull
10-0570 for patient surface pull-out screwed tight
10-0628 for patient surface width of 100 cm and patient surface pull-out screwed tight



10-0563/10-0533/10-0521
Segufix holder for thigh rest
Segufix holder for the backrest
Segufix holder for the seat section



10-0595
Patient surface extension for beds with
integrated extension, wooden decor



**Bed extension
for beds with integrated patient surface
extension**
10-0497
Extension element to extend the lower
leg rest by 20 cm
10-0573
Extension element to extend the calf rest
by 20 cm with 100 cm wide patient surface

Protection

Depending on how the bed is used, it is possible to prevent the resident from injuries and the bed from damage caused by the surroundings. Side guard cushioning compliments the bed and closes the gaps.



50-0206
Fall mat optimal for sentida beds



10-0597
Side guard extension
for continuous (09) side guard

10-0605
Side guard extension
for divided (10) side guards (no image)



50-0544
Evacuation sheet with belt and fixation straps



50-0155
Side section padding with push buttons
for divided (10) side guards



50-0174
Side section padding for divided (10)
side guards with safety strap, closed type



50-0207
Side guard net to close the gaps
of divided (10) side guards



50-0338
Side section padding for continuous (09)
two-part side guards for bed length of 200 cm
50-0339
Side section padding for continuous (09)
two-part side guards for bed length of 220 cm
(no image)

Protection



10-0682
3D wall stopper castor, reinforced type
10-0683
3D wall stopper castor for design A (no image)



10-0485
Reinforced wall stopper castors for divided (10) side guards as well as continuous 2-part and 3-part 09 side guards




10-0626
Vertical wall stopper castors



Vertical wall stopper castors, extended type
10-0507 extended by 22 cm
19-0494 extended by 12 cm
19-0480 extended by 9 cm



02-000617
Side panel buffer for divided (10) side guards



Side-mounted wall spacer
10-0549 for movita/movita casa
10-0578 for carisma



10-0530
Head-end wall spacer bracket carisma
Distance to wall approx. 17.4 cm



10-0548
Head-end wall spacer bracket movita
Distance to wall approx. 8 cm



10-0587
Head-end wall spacer bracket movita casa
Distance to wall approx. 11.5 cm



Head-end wall spacer, telescopic, only for sentida 5-le and 6-le
10-0632 short, 1 x 20 mm
10-0622 medium, 4 x 20 mm
10-0612 long, 8 x 20 mm



Side-mounted wall spacer, telescopic, only for sentida 5-le and 6-le
10-0633 short, 1 x 20 mm
10-0623 medium, 4 x 20 mm
10-0613 long, 8 x 20 mm



Head-end wall spacer, telescopic, only for sentida 5 and 6
10-0801 short, 1 x 20 mm
10-0802 medium, 4 x 20 mm
10-0803 long, 8 x 20 mm



Side-mounted wall spacer, telescopic, only for sentida 5 and 6
10-0804 short, 1 x 20 mm
10-0805 medium, 4 x 20 mm
10-0806 long, 8 x 20 mm

Head-end wall spacer, only for sentida 5 and 6 (no image)
10-0707 short
10-0709 long

Side-mounted wall spacer, only for sentida 5 and 6 (no image)
10-0708 short
10-0710 long

Care

The accessories for the accessory holder provide useful support for patient and staff.
The holder can be attached to the left and right side of the bed.



50-0693
Bent infusion stand with 4 hooks,
height-adjustable



50-0003
Drip clasp to attach to the lifting pole



50-0694
Straight infusion stand with 4 hooks
and drip vessel **01-006094**



10-0534
Transport handle „Easy Move“



50-0108
Name tag holder



10-0424
Accessories holder for divided (10) side
guards, to attach to the side panel



50-0345
Shortened urine bottle basket
similar to picture)



50-0164
Urine bag holder (depicted here with
urine bag and accessories holder)
59-0031
Extended urine bag holder (no image)

10-0511
Bed transport trolley for assembled
nursing home beds
10-0536
Bed transport trolley, type B,
for unassembled care beds (no image)

Members of LINET Group

wissner-bosserhoff GmbH

Hauptstraße 4-6 | 58739 Wickede (Ruhr) | Germany

Phone: +49 2377.784-0 | Fax: -163 | info@wi-bo.de | wi-bo.com

