- wissner-
- bosserhoff



Elderly care

Maximum functionality meets modern design























wissner-bosserhoff products offer more than the average nursing home bed, bedside cabinet or mattress. In cooperation with nursing staff, added value is identified, developed and offered to the user.

wissner-bosserhoff stands for:

- Digital and networked nursing home beds
- Designer upholstery beds with hotel-like character
- Unique and proven solutions for mobilization, pressure ulcer and falls prevention

Ever since we were founded in 1952, we have been a family-owned company and have maintained long-standing partnerships with our customers.

Convince yourself!

Your wissner-bosserhoff - team

360° Elderly care 2 | 3









Many older people find it difficult to leave their homes when sufficient health and nursing care can no longer be ensured at home. This makes it even more challenging for nursing homes to offer these residents a new and comfortable home.

In contrast to most nursing home beds with a predominantly functional design, wissner-bosserhoff excels through its outstanding homeliness, and thus actively contributes to a "healing environment".

Our solutions

- Homely ambience for an active healing environment
- New upholstered bed collection in an attractive hotel design, with a vast number of fabric and imitation leather designs
- New, homely side guard concept SafeFree® Flex
- Extensive range of design and decor options available
- Numerous award-winning designs and products in recent years

More about our upholstered beds in hotel design:














360° Elderly care 4 | 5



Nursing crisis is already having a real and palpable impact. In many countries the society is aging rapidly and the number of care-dependent people is set to increase rapidly in the coming years. In contrast, more and more nursing facilities are in need of skilled workers.

New digital options not only relieve the burden on nursing care by significantly saving time, they also contribute to greater cost-effectiveness.

Our solutions

- Direct connection of our sentida 7-i nursing home bed to the nursing documentation
- All major interface standards integrated Secure and reliable real-time data transfer
- Calibrated, integrated weighing scale in sentida 7-i with integrated BMI check
- Differentiated nurse call evaluation for correct prioritization
- sentida 7-i early warning system for the fastest possible response to detectable fall cases

More about the sentida 7-i digital nursing home bed:



71,8 kg 71 5 6 7 7 1 8 kg 7 8 1 8 kg 7 8 kg

Nursing documentation



Integrated scale



Differentiated nurse calls



Bed-Exit-Alarm


360° Elderly care 6 | 7

Falls prevention









Falls pose an increasing, overall risk as residents grow older. Statistically speaking, 30 % of all people above 65 fall at least once a year. In the case of people over 80, this figure increases to almost 50 %. The percentage of people at risk of falling is even higher in nursing homes. According to statistics, more than half of residents fall at least once per year.

Source: www.sturzpraevention.net

Our solutions


- SafeSense® bed-exit assistant for effective falls management with automatic night light for safe orientation in the dark
- Flexible nurse call timing for faster response time in the event of falls
- SafeFree® side guards, individually adjustable according to care situation
- Low position on all wissner-bosserhoff care beds for falls prevention during sleeping phases
- Comprehensive safety sensor technology, e.g. brake assist or side guard assist in the nursing home bed sentida 7-i


More about the SafeSense® Bed Exit - System:











SafeSense®

Night light

SafeFree®

Safety sensors



In nursing care, maintaining the mobility of residents in the long term is a major challenge. Movement and mobility have a positive and preventive effect on many areas of everyday care. However, the mobilization of residents in nursing care facilities does not always come without the risk of falls. Our products can support residents and provide them with additional safety.

Our solutions

- Revolutionary, wireless mobilization assistant SafeLift with simplest operation (height adjustment and backrest) for the resident*
- MobiLift® mobilisation handle in the sentida 7-i nursing home bed with integrated height adjustment
- New, manual MobiStick2 mobilization-assistant with intuitive operation*
- Split SafeFree® side guards as a safe standing aid, individually adaptable to different heights of residents
- 3-stop strategy with automatic intermediate stop at ideal mobilization height

More about the SafeLift Mobilization Assistant:











SafeLift

MobiLift®

MobiStick2

SafeFree®

^{*} for all wissner-bosserhoff care beds with SafeFree® side guards



Relocating, moving and mobilizing patients are integral parts of a nurse's daily work. Anyone working with the elderly and patients knows: This is back-breaking work in the truest sense of the word. Many nursing staff therefore suffer from back problems. This is most often caused by frequent lifting of heavy loads - often in a forward bent posture - and lack of use of technical aids in the fast-paced and often stressful working day.

Our solutions

- Optimal nursing position of 80 cm for back-friendly execution of care routines
- Clear and ergonomic controls
- Safety through compliance with the IEC 60601-2-52 standard for the safety of medical beds.
- Integrated bed extension as standard in all nursing beds.
- carisma 300-xxl nursing home bed for the safe care of obese residents

More about the nusing bed for obese residents carisma 300-xxl:



Nursing position



Clear controls



Integrated bed-extension



carisma 300-xxl

360° Elderly care 12 | 13



Pressure ulcers are very unpleasant for residents and can affect their quality of life. They suffer from severe pain and treatment usually lasts several weeks or even months. Care facilities may incur considerable follow-up costs. That's why the development of pressure ulcers should be avoided altogether wherever possible. Preventive measures can be taken through the use of aids.

Our solutions

- Unique BiNetic patient surface with ergonomic layout
- Double retraction as part of the patient surface to reduce pressure loads
- Comfort patient surface with 55 spring elements
- Extensive range of passive mattresses, for optimum pressure distribution at all lying body parts
- Large selection of active and reactive therapy systems for pressure ulcer prevention and therapy with state-of-the-art technologies

More about our extensive range of mattresses:




0

WIBO patient surface with notable relief



Double retraction = pressure ulcer prevention



Thermal image of a BiNetic patient surface



Extensive range of mattresses

360° Elderly care 14 | 15



Member of LINET Group

wissner-bosserhoff GmbH

Hauptstr. 4-6 | 58739 Wickede (Ruhr) | Germany

Phone: +49 2377 784-0 | Fax: -167 | info@wi-bo.de | www.wi-bo.de