




SafeLift Intuitive simplicity





Improving and maintaining mobility

SafeLift | Revolutionary bed operation for residents

A key goal in in-patient care for the elderly is to maintain or improve mobility. More mobile residents are less prone to falling, further secondary diseases such as pressure ulcers are prevented, and their health is kept in optimum condition. This increases their self-esteem and an independent lifestyle is less restricted.

SafeLift offers the simplest operation so that residents can adjust their bed independently and set it to the best mobilization position for them. At the same time, they can support themselves when standing up on the SafeLift controller.



Simple bed operation for residents



Standing-up and mobilization aid



Cable- and battery-free controller



Customized for each resident

wissner-bosserhoff | SafeLift

Simple bed operation & safe mobilization

SafeLift | Product benefits at a glance





Adjustment buttons for bed height and backrest only the features which residents really needs



Large buttons easy to feel and fluorescent lighting



Firm position on a SafeFree® side guard make it easier to find



GO button protects against unintentional bed adjustment



360° rotation and flexible adjustment on the SafeFree® side guard



Ergonomic design suitable for any size of hand and easy accessibility of the operating buttons



Strong stability provides firm grip when standing up





Residents can set their bed independently to the best position to get out and support themselves on the SafeLift controller.



No cable or battery

Thanks to established radio technology, the SafeLift controller can be operated without any fault-prone cables.



Automatic night light

and fluorescent buttons make it easier to get around in the dark.









4 wissner-bosserhoff | SafeLift wissner-bosserhoff | SafeLift

SafeLift – intuitive simplicity

Added value for residents, care and care-home management



Resident





- Accessible to residents when lying down, individually positionable and usable as a support
- Fall prevention no tripping hazards from hanging cables
- Position of the SafeLift controller individually adjustable for each resident
- Suitable for right- and left-handed people, as well as people who get out of the right or left side of the bed
- Automatic night light makes it easier to get around in the dark



Nursing staff



Promotes relief & mobility in line with expert standards

- Promotes the mobility and independence of the resident
- For routine tasks, always ready to hand for care, too, care positions as usual by hand switch
- Automatic night light
- Simple and tool-free installation on every SafeFree® side guard*
- Extremely easy to adjust to individual resident needs



Management



Flexibility and economic efficiency

- Low follow-up costs –
 no fault-prone cables and no battery
- SafeLift controller can be attached without tools to any SafeFree® side guard
- Can be used on any bed
- Compatible with the SafeSense®

 Bed Exit Assistance System
- Can be retrofitted to all wissner-bosserhoff care beds with SafeFree® side guards*







* Exceptions: carisma 300, carisma 300-xxl, sentida 1p/4p, sentida 7i wissner-bosserhoff | SafeLift wissner-bosserhoff | SafeLift





Members of LINET Group

wissner-bosserhoff GmbH

Hauptstr. 4-6 | 58739 Wickede (Ruhr) | Germany

Phone: +49 2377.784-0 | Fax: -163 | info@wi-bo.de | wi-bo.com

